

QPCR Reagent Selection Guides

ABsolute Fast QPCR - ABsolute Blue QPCR - ABsolute QPCR - Verso QRT-PCR

ABsolute Fast QPCR Master Mix Selection Guide

QPCR Cycler	ABsolute Fast QPCR Master Mix Plus ROX Vial	ABsolute Fast QPCR ROX Master Mix	ABsolute Fast QPCR Low ROX Master Mix
ABI PRISM®			
7000/7300		•	
7500			•
7700		•	
7900/7900HT		•	
StepOne™/StepOnePlus™			•
Bio-Rad			
iCycler™	•		
MyiQ™	•		
iQ™ 5	•		
Opticon™/2	•		
Chromo 4™	•		
MiniOpticon	•		
CFX 96	•		
Stratagene			
Mx4000™	•		•
Mx3000™/3005™	•		•
Techne			
Quantica™	•		
Cepheid			
SmartCycler™	•		
Corbett			
Rotor-Gene™	•		
Roche			
Lightcycler® 480	•		
Eppendorf			
Realplex	•		

ABsolute Blue QPCR Master Mix Selection Guide

QPCR Cycler	P R O B E			S Y B R			
	Mix + ROX Vial	ROX Mix	Low ROX Mix	SYBR Green Mix + ROX Vial	SYBR Green ROX Mix	SYBR Green Low ROX Mix	SYBR Green Fluorescein Mix
ABI							
7000/7300		•			•		
7500			•			•	
7700		•			•		
7900/7900HT		•			•		
StepOne™/StepOne™ Plus			•			•	
Bio-Rad							
iCycler™	•						•
MyiQ™	•						•
iQ™ 5	•			•			
Opticon™/Opticon™ 2	•			•			
Chromo 4™	•			•			
MiniOpticon	•			•			
Stratagene							
Mx4000™	•		•	•		•	
Mx3000P™/3005P™	•		•	•		•	
Techne							
Quantica™	•			•			
Cepheid							
SmartCycler™	•			•			
Corbett							
Rotor-Gene™	•			•			
Roche							
Lightcycler® 480	•			•			
Eppendorf							
Realplex	•			•			

Absolute QPCR Master Mix Selection Guide

QPCR Cycler	P R O B E					S Y B R					
	Mix + ROX Vial	ROX Mix	Low ROX Mix	Mix (No ROX)	Capillary Mix (No ROX)	SYBR Green Mix + ROX Vial	SYBR Green ROX Mix	SYBR Green Low ROX Mix	SYBR Green Mix (No ROX)	SYBR Green Fluorescein Mix	SYBR Green Capillary Mix (No ROX)
ABI											
7000/7300		•					•				
7500			•					•			
7700		•					•				
7900/7900HT		•					•				
StepOne™/StepOne™ Plus			•					•			
Bio-Rad											
iCycler™	•			•						•	
MyiQ™	•			•						•	
iQ™ 5	•			•		•			•		
Opticon™/Opticon™ 2	•			•		•			•		
Chromo 4™	•			•		•			•		
MiniOpticon	•			•		•			•		
Stratagene											
Mx4000™	•		•			•		•			
Mx3000P™/3005P™	•		•			•		•			
Techne											
Quantica™	•			•		•			•		
Cepheid											
SmartCycler™	•			•		•			•		
Corbett											
Rotor-Gene™	•			•		•			•		
Roche											
Lightcycler® 2.0					•						•
Lightcycler® 480	•			•		•			•		
Eppendorf											
Realplex	•			•		•			•		

Verso QRT-PCR Master Mix Selection Guide

QPCR Cycler	P R O B E			S Y B R			
	+ ROX Vial	Incl. ROX	Low ROX	SYBR Green + ROX Vial	SYBR Green ROX	SYBR Green Low ROX	SYBR Green Fluorescein
ABI							
7000/7300		•			•		
7500			•			•	
7700		•			•		
7900/7900HT		•			•		
StepOne™/StepOne™ Plus			•			•	
Bio-Rad							
iCycler™	•						•
MyiQ™	•						•
iQ™ 5	•			•			
Opticon™/Opticon™ 2	•			•			
Chromo 4™	•			•			
MiniOpticon	•			•			
Stratagene							
Mx4000™	•		•	•		•	
Mx3000P™/3005P™	•		•	•		•	
Techne							
Quantica™	•			•			
Cepheid							
SmartCycler™	•			•			
Corbett							
Rotor-Gene™	•			•			
Roche							
Lightcycler® 480	•			•			
Eppendorf							
Realplex	•			•			